




















## Non-Implantable Bone Conduction Devices

		Suitable for babies 	Suitable for Unilateral 	Suitable for Bilateral 	How attached 	Can device be implanted later? 
Oticon Medical Ponto		Yes	Yes	Yes	Headband	Yes
Cochlear Baha		Yes	Yes	Yes	Headband or SoundArc band	Yes
MED-EL ADHEAR		Yes	Yes	Yes	Adhesive adapter or Headband	No

## Implantable Devices

		How device works 	How processor is attached 	Suitable for Unilateral / Bilateral 	Waterproof option 	Approved age for implant 
Oticon Medical Ponto		Bone conduction	Abutment	Yes	No	5 years
Cochlear Baha Attract		Bone conduction	Magnet	Yes	No	5 years
Cochlear Baha Connect		Bone conduction	Abutment	Yes	No	5 years
Cochlear Osia		Active Implant	Magnet	Yes	Yes	10 years
MED-EL BONEBRIDGE		Active Implant	Magnet	Yes	Yes	5 years
MED-EL SOUNDBRIDGE		Active Implant	Magnet	Yes	Yes	5 years